

Version	New functions	Known bugs	What happens?	Remarks / Infos	Date
Netus 3.2.12 / TraceMon 3.2.12	New method "Dialysis Low" is available	Measurement external sample in MultiTRACE mode	The reference factor will not be assigned to the right probe	Program from 16.01.2017	16.01.2017
Netus 3.2.15		Known bugs from former versions solved		Updater from 05.11.2018	05.11.2018
Netus 5.1.20 / TraceMon 5.1.20	Data transfer via OPC-UA available. Additional device M4 can be connected.	Known bugs from former versions solved		Release-Updater from 13.12.2019	13.12.2019
TraceMon 5.1.21	Russian Language optimized			Program from 07.10.2020	07.10.2020
TraceMon 5.1.22		Bug in Autoreference for C2/M4 Methanol/Ethanol fixed		Program from 20.04.2021	20.04.2021